

## All in One™ Garment Sizing Chart

Sleeve	Size	Catalog Number	Proximal Circumference	Distal Circumference
Calf Length (single chamber)	S	L4ZD261M	29-45	20-30
Calf Length (single chamber)	M	L4ZD262M	34-60	24-36
Calf Length (single chamber)	L	L4ZD263M	46-64	34-47
Calf Length (four chamber)	S	L4ZD211M	29-45	20-30
Calf Length (four chamber)	M	L4ZD212M	34-60	24-36
Calf Length (four chamber)	L	L4ZD213M	46-64	34-47
Calf and Foot (four chamber)	S	L4ZD2310	32-42	27-33
Calf and Foot (four chamber)	M	L4ZD2320	36-46	26-36
Calf and Foot (four chamber)	L	L4ZD2330	42-54	36-44
Thigh Length (four chamber)	S	L4ZD2410	44-60	24-28
Thigh Length (four chamber)	M	L4ZD2420	54-76	32-42
Thigh Length (four chamber)	L	L4ZD2430	54-82	30-47
Foot Cuff	M	L4ZD272M	UP TO SIZE 39	

Manufactured by:  
Mego Afek | Kibbutz Afek 3004200 Israel  
Tel: +972-77-9084277 Fax: +972-4-8784148

[info@meoafek.com](mailto:info@meoafek.com)

## Phlebo Press® DVT 603

*All in One™*

The Hospital Choice for DVT Prevention



Designed to meet clinician and patient needs.

Phlebo Press® DVT All in One™ is fully automatic, for easy operation and good compliance.

### Excellence in DVT Prevention

- Easy to use – On/Off operation! Automatically operates any sleeve type and combination, at correct pressure.
- Cost effective – Use with single or four chamber sleeves, calf, thigh or foot.
- Battery power – Up to eight hours of battery operation.
- Stable bed rail mount – For easy nurse access.
- Exceptionally soft, comfortable sleeves





Easy to use!  
Attach any sleeve type,  
press START, and the  
system does the rest.

Choice of gradient circumferential compression with four-chamber sleeves, or economical asymmetric single-chamber compression in one device.

- Gradient sequential circumferential compression at 45 mmHg with four chamber sleeves.
- Intermittent compression with single chamber sleeves.
- Foot sleeves operate at 130 mmHg – for orthopedic procedures requiring access to knee area.
- Automatic recognition for any sleeve combination – one or both legs, foot, calf or thigh. Press START and the device identifies and applies correct pressure and therapy.
- Battery power for up to eight hours of mains-free use, recharges continuously during mains operation.
- Clear audible and visual alarms. Easy identification of common use errors including kinked hose, disconnected sleeve or improper garment fastening.
- Sleeves are soft and easy to wear. Unique ventilated sleeve materials keep the leg cool and comfortable. Sleeves direct hose away from ankle and heel, preventing irritation and skin breakdown.

The *All in One*<sup>™</sup> identifies sleeve type  
– any combination - automatically.  
Side lights make it easy to monitor  
operation.

	POWER BUTTON		BATTERY STATUS
	START/STOP BUTTON	130	PRESSURE LEVEL (mmHg)
	GARMENT DISCONNECT/AIR LEAK		CONNECTED GARMENT(S)
	KINKED HOSE		STATUS INDICATOR Green: Treatment in progress. White: Device powered, treatment not running. Red: Device in error state.
	AC POWER		
	BILATERAL TREATMENT	4	SINGLE OR QUAD CHAMBERS
CYCLE 60 seconds	WEIGHT AND SIZE 2.3 kg / 23 x 20 x 18 in.	PRESSURE RANGE 45 / 130 mmHg	ELECTRICAL 90-246 VAC, 47-63 Hz

## DESIGNED FOR COMPLIANCE

### Automatic operation

The *All in One*<sup>™</sup> automatically detects and operates any sleeve combination.

### Easy to use and monitor

Sleeve shape and markings help ensure correct application on every patient. Side lights allow the caregiver to see operating status at a glance.

### Soft, wearable sleeves

Soft, ventilated sleeves keep patient's legs comfortable and cool and are designed to prevent pressure points and skin breakdown.

